

DIGITAL POOL & SPA TESTING

# POOL CHECK™

## eXact™ Strip

**PART NUMBER 481647**

The eXact™ Strip Cl Free (481647) offers the safest and most convenient delivery system for your Free Chlorine testing with the eXact™ Meter. All reagents, as identified in the Standard Methods AWWA Handbook (4500-Cl G.), are contained in the eXact™ strip pads.



### **Pool Check™ eXact™ Strip For Free Chlorine**

**(For Use With eXact™  
1,2,3,& 9 Meters)**

Available For Private Label.

#### **Test Procedure (eXact™ 1 meter example):**

**Step 1:** Remove one eXact™ Strip Cl Free from the bottle. Recap the bottle immediately. Set strip in a dry place.

**Step 2:** Rinse supplied plastic vial 1 or 2 times and fill to line with pool water sample. DO NOT OVERFILL! Dry outside of vial with a soft paper or cloth towel.

**Step 3:** Lift the chamber lid and insert plastic vial. Close lid.

**Step 4:** Press the "Menu" key to switch the meter on. Press the "Menu" key to display the "Cl" symbol.

**Step 5:** Press the "B" key to blank the meter to 0.

**Step 6:** Open the chamber lid and dip the eXact™ Strip Cl Free into the plastic vial for 20 seconds with gentle back and forth motion.

**Step 7:** Remove and discard the strip; immediately close the lid and push the "R" key. Record result.

#### **Detection Range:**

**Free Chlorine: 0-5.00 mg/L**

#### **Packaging:**

**Bottle of 50 eXact™ strips**

## **INDUSTRIAL TEST SYSTEMS, Inc.**

**1875 Langston Street, Rock Hill, SC 29730**

**Phone: 1-800-861-9712 Fax: 1-803-329-9743**

**For International: 001-803-329-9712 ext 220 (Pablo)**

**Web: [www.poolcheckonline.com](http://www.poolcheckonline.com) Email: [its@cetlink.net](mailto:its@cetlink.net)**



Member of



**APSP**

The Association of  
Pool & Spa Professionals